

VOMIT BAG

The vomit bag offers a convenient and simple way to collect and dispose of vomit

- Wide, convenient intake opening to prevent spillage of contents;
- Made entirely of transparent material that allows you to visually examine the contents inside the bag, including blood
- The plastic rim allows easy grip and handling
- Rugged bag with a capacity of 1.5L
- with scale (ml)
- The „Twist & lock” function on the rim allows you to securely close the bag after use
- Practical, easy to use
- Available version with liquid absorber

REF	FLUID ABSORBER	PACKAGING
06001	NO	50 pcs.
06001Z	YES	50 pcs.

